

Jigsaw Medical®

Endoscopy Product Catalog

Procedure Specific Solutions — *Carpal Tunnel Release*

Jigsaw Medical® ECTR Blade Kit

Our sterile, single-use kit is compatible with Smith+Nephew™ ECTRA® II System to deliver reproducible results in the precise access and release of the carpal ligament while simultaneously improving the bottom line.

Contents:

- 1 Triangle Blade (green)
- 1 Retrograde Blade (blue)
- 1 Probe Blade (natural)
- 10 Cotton Applicators



Ref #	Description	UOM
50119-S	Sterile, Single-Use ECTR Kit	each



**Jigsaw
Medical®**

Jigsaw Medical, LLC
1012 Airpark Dr. Suite 1
Sugar Grove, IL 60554



www.jigsawmed.com



(630) 587-0019



info@jigsawmed.com

Jigsaw Medical® Individual ECTR Blades

Triangle Blade with Green Handle



- Sterile, single-use individual blade
- Incises the midsection of the carpal ligament to create an opening for insertion of the Retrograde Blade.
- Sharp tip of the Triangle Blade allows for exactness in follow-up cuts dissecting individual fibers.

Ref #	Description	UOM
50116-S	Triangle Blade with Green Handle	6/box

Retrograde Blade with Blue Handle



- Sterile, single-use individual blade
- Precisely cuts the carpal ligament without damaging the palmaris brevis.
- Sharp upper-tip of the Retrograde Blade smoothly glides between ligament and muscle while the ligament is guided into the blade's lower-arch.

Ref #	Description	UOM
50118-S	Retrograde Blade w/ Blue Handle	6/box

Probe Blade with Natural Handle



- Sterile, single-use individual blade
- Releases the distal aspect of the ligament while providing protection to surrounding neurovascular structures.
- Rounded, blunt top of the Probe Blade dissects synovium to precisely locate ligament's edge.

Ref #	Description	UOM
50117-S	Probe Blade with Natural Handle	6/box



Jigsaw Medical, LLC
1012 Airpark Dr. Suite 1
Sugar Grove, IL 60554

 www.jigsawmed.com
 (630) 587-0019
 info@jigsawmed.com