



Product Catalog

Pleural and Peritoneal Fluid Evacuation



Jigsaw Medical™

OVERVIEW

Jigsaw Medical® Drainage Products

Jigsaw Medical® carries the broadest range of evacuated drainage bottles on the market. Our Safety-Vac™ line consists of four evacuated drainage bottles and two drain lines for optimal performance in removing, aspirating, and sampling body fluid and assured compatibility with other market-leading centesis drainage products. Jigsaw Medical®'s specialization in this area of healthcare products has rapidly become a clinically-preferred solution to providing high quality care, while simultaneously improving the bottom line.

1000+mL Drainage Bottles

Spike-able Sterile Pathway Safety-Vac™
(SV1000-SSP)

Standard Sterile Safety-Vac™ + Tubing
(SV1000-S)

Standard Sterile Pathway Safety-Vac™
(SV1000-NS)

Valved Standard Sterile Pathway Safety-Vac™
(SV1000-NSWC)

Drainage Accessories

High-Flow Sterile 47" inch Drain Line
(SVHFSP-S)

Standard Sterile 47" inch Drain Line
(SVCT47-S)



NOT MADE WITH NATURAL
RUBBER LATEX



NOT MADE WITH DEHP



NOT MADE WITH PVC

Spike-able Sterile Pathway Safety-Vac™

1000+mL Evacuated Drainage Bottle

PATENT PENDING

- Widely Used for Paracentesis, Thoracentesis, and Other Drainage Procedures.
- Intended for Use to Aspirate, Remove, or Sample Body Fluid.
- Not Made with Natural Rubber Latex, DEHP, or PVC.
- Vacuum up to 27inHg (Inches of Mercury) of Negative Pressure.

Feature	Benefit
Dual-Functioning	Compatible with needle-tip AND luer-lock drain lines
1325mL Max. Fill	Higher than glass bottles & standard plastic bottles
Made of high-quality, medical grade plastic	Shatter-resistant and highly transparent for visibility
Sterile Pathway <i>(internally sterilized)</i>	Accurate lab samples
12 Bottles/Case	Increased case quantity equals cost savings
Exclusive only to Jigsaw Medical®	No other company supplies this product



REF #



SV1000-SSP

ITEM

Spike-able Sterile Pathway Safety-Vac™
1000+mL Evacuated Drainage Bottle

CASE QTY.

12



Standard Sterile Safety-Vac™ + Tubing

1000+mL Individually Packaged Evacuated Drainage Bottle




PATENT PENDING

- Widely Used for Paracentesis, Thoracentesis, and Other Drainage Procedures.
- Intended for Use to Aspirate, Remove, or Sample Body Fluid.
- Sterile Field Ready
- Not Made with Natural Rubber Latex, DEHP, or PVC.
- Vacuum up to 27inHg (Inches of Mercury) of Negative Pressure.

Feature	Benefit
Internal/External Sterilization	Sterile field ready; Accurate lab samples
Sterile Luer-Lock Tubing Included	Needle-free drainage; Additional cost-savings
Made of high-quality, medical grade plastic	Shatter-resistant and highly transparent for visibility
Holds 1000+mL of Fluid	Slightly higher than market standards
12 Bottles/Case	Increased case quantity equals cost savings



REF #	ITEM	CASE QTY.
SV1000-S	Standard Sterile Safety-Vac™ + Tubing 1000+mL Evacuated Drainage Bottle	12



Standard Sterile Pathway Safety-Vac™

1000+mL Evacuated Drainage Bottle


PATENT PENDING

- Widely Used for Paracentesis, Thoracentesis, and Other Drainage Procedures.
- Intended for Use to Aspirate, Remove, or Sample Body Fluid.
- Not Made with Natural Rubber Latex, DEHP, or PVC.
- Vacuum up to 27inHg (Inches of Mercury) of Negative Pressure.

Feature	Benefit
Luer-Lock Compatible	Needle-Free Drainage
Holds 1000+mL of Fluid	Slightly higher than market standards
Made of high-quality, medical grade plastic	Shatter-resistant and highly transparent for visibility
Sterile Pathway <i>(internally sterilized)</i>	Accurate lab samples
12 Bottles/Case	Increased case quantity equals cost savings



REF #	ITEM	CASE QTY.
SV1000-NS	Standard Sterile Pathway Safety-Vac™ 1000+mL Evacuated Drainage Bottle	12



NOT MADE WITH NATURAL RUBBER LATEX



NOT MADE WITH DEHP



NOT MADE WITH PVC



SAFETY-VAC™ BOTTLES

Valved Standard Sterile Pathway Safety-Vac™

1000+mL Evacuated Drainage Bottle




PATENT PENDING

- Widely Used for Paracentesis, Thoracentesis, and Other Drainage Procedures.
- Intended for Use to Aspirate, Remove, or Sample Body Fluid.
- Not Made with Natural Rubber Latex, DEHP, or PVC.
- Vacuum up to 27inHg (Inches of Mercury) of Negative Pressure.

Feature	Benefit
Valved Luer-Lock Connection	Needle-free, completely closed suction system
Holds 1000+mL of Fluid	Slightly higher than market standards
Made of high-quality, medical grade plastic	Shatter-resistant and highly transparent for visibility
Sterile Pathway <i>(internally sterilized)</i>	Accurate lab samples
12 Bottles/Case	Increased case quantity equals cost savings



REF #	ITEM	CASE QTY.
SV1000-NSWC	Valved Standard Sterile Pathway Safety-Vac™ 1000+mL Evacuated Drainage Bottle	12




Jigsaw Medical® Drainage Accessories

- Intended for Use to Aspirate, Remove, or Sample Body Fluid Using an External Vacuum Source.
- Not Made with Natural Rubber Latex, DEHP, or PVC.

High-Flow Sterile Drain Line

SVHFSP-S





Feature	Benefit
Wider, thicker bore	30% flow speed increase using external vacuum sources
Optional Plastic Spike Attachment	Can be used with glass bottles
Fast Drainage Speed	Less time spent in high volume paracentesis drains
Sterile, 47 inch Luer-Lock Drain Line	Optimal compatibility with either a bottle or wall suction
50 Drain Lines/Case	Typically at a lower cost than market standards



Standard Sterile Drain Line

SVCT47-S

Feature	Benefit
Attached Roller-Clamp	Easily controllable flow rate from external vacuum sources
Luer-Lock Compatibility	Needle-free drainage
Sterile, 47 inch drain line	Optimal compatibility with either a bottle or wall suction
50 Drain Lines/Case	Typically at a lower cost than market standards



ORDERING INFORMATION

Safety-Vac™ Evacuated Drainage Bottles

REF #	ITEM	CASE QTY
SV1000-SSP	Spike-able Sterile Pathway 1000mL Evacuated Drainage Bottle	12
SV1000-S	Standard Sterile 1000mL Evacuated Drainage Bottle + Tubing	12
SV1000-NS	Standard Sterile Pathway 1000mL Evacuated Drainage Bottle	12
SV1000-NSWC	Valved Standard Sterile Path. 1000mL Evacuated Drainage Bottle	12

Safety-Vac™ Evacuated Drainage Bottles are intended for use to aspirate, remove, or sample body fluids.

Drainage Accessories

REF #	ITEM	CASE QTY
SVHFSP-S	High-Flow Sterile 47" inch. Drain Line w/ Optional Spike Attachment	50
SVCT47-S	Standard Sterile 47" inch. male-male Drain Line	50



Jigsaw Medical® Drainage Accessories intended use: To aspirate, remove, or sample body fluids using an external vacuum source.

For all order inquiries, pricing, or product questions, please contact your Jigsaw Medical® representative.

Tel (630) 587-0019 **Fax** (630) 524-9088

Email info@jigsawmed.com

Refer to Instructions for Use for relevant warnings, precautions, complications, and contraindications.



Jigsaw Medical®
1012 Airpark Dr. Suite 1
Sugar Grove, IL 60554
www.jigsawmed.com
(630) 587-0019