


Exclusive only to Jigsaw Medical®

# TrapDoor™ Wound Cover

Effective Solution in Wound Care Management



PATENTED

## The Only Re-Closeable Wound Cover On The Market.

Through innovative design and prioritizing patient care, the Jigsaw Medical® TrapDoor™ Wound Cover is a revolutionary solution to efficient daily wound management for patients and providers. The re-closeable TrapDoor Wound Cover can be worn for 7 continuous days without removing any adhesive from the patient's skin, which greatly improves the overall wound healing process.


### Wear for 7 Days

The breathable primary base adhesive can be worn for 7 continuous days. Longevity benefits also include:

- **Reduced Frequency of Painful Film/Bandage Removal**
- **Reduced Skin Tear Risk**
- **Faster Wound Healing**
- **Improved Patient Comfort**
- **Time + Cost Efficient**

### Size Variety

TrapDoor™ Wound Cover comes in five (5) different sizes to appropriately protect all wounds.



### Water-Resistant

TrapDoor™ Wound Cover is water-resistant and can be worn in the shower if direct flow is avoided.

### Re-Closeable Window

Allows for daily access to the wound/port site without removing any adhesive from the skin. Other benefits include:

- **Durable, Flexible Window Flap**
- **Convenient Daily Wound Access**
- **Protects + Stabilizes Dressing**
- **Easy Catheter Port Access**
- **Reduced Frequency of Painful Film/Bandage Removal**

**7 Days.**  
**5 Sizes.**  
**Zero Removal.**

TrapDoor™ Wound Cover is  
**Revolutionizing**  
**Wound Care.**

REF #	ITEM	UOM
TD0605	X-Small TrapDoor Wound Cover	25/box
TD1310	Small TrapDoor Wound Cover	10/box
TD1822	Medium TrapDoor Wound Cover	10/box
TD3012	Large TrapDoor Wound Cover	10/box
TD3023	X-Large TrapDoor Wound Cover	10/box



## TrapDoor™ Size Chart

Not intended for use as a primary wound closure device.  
Product pictures are not to scale.



### TD0605 | X-Small


Dressing Size: 2.5" x 2.2"  
Window Size: 1" x 1"  
Box Quantity: 25

*Box includes dressing change chart.*

### TD1310 | Small

Dressing Size: 5.1" x 4"  
Window Size: 2.25" x 2.25"  
Box Quantity: 10

*Box includes dressing change chart.*



### TD1822 | Medium


Dressing Size: 7" x 8.6"  
Window Size: 4" x 4"  
Box Quantity: 10

*Box includes dressing change chart.*

### TD3012 | Large

Dressing Size: 11.8" x 4.7"  
Window Size: 8" x 2.3"  
Box Quantity: 10

*Box includes dressing change chart.*



### TD3023 | X-Large

Dressing Size: 11.8" x 9"  
Window Size: 6" x 6"  
Box Quantity: 10

*Box includes dressing change chart.*

For all order inquiries, pricing, or product questions, please contact your Jigsaw Medical® representative.

**Tel** (630) 587-0019    **Fax** (630) 524-9088

**Email** [info@jigsawmed.com](mailto:info@jigsawmed.com)

Refer to Instructions for Use for relevant warnings, precautions, complications, and contraindications.



**Jigsaw Medical®**  
1012 Airpark Dr. Suite 1  
Sugar Grove, IL 60554  
[www.jigsawmed.com](http://www.jigsawmed.com)  
(630) 587-0019