

Exclusive only to Jigsaw Medical®

The Spike-able Safety Vac™

2-in-1 Drainage Bottle



*Patent Pending

SV1000-SSP

The Only Plastic Spike-able Evacuation Bottle on the Market.

Through advanced molding capabilities, the Jigsaw Medical® **Spike-able Safety Vac™** combines the functionality of glass and plastic evacuated bottles to meet all clinical preferences while also improving overall hospital safety, patient care quality, and the bottom line.

Dual-Functionality

Clinicians can choose their preferred drainage method:

Method 1: Spike cork with a needle tipped drain line.

Method 2: Connect luer lock drain line.

Sterile Pathway

Internally sterilized for accurate lab results.

Easy Sample Access

Insert syringe through cork to easily retrieve samples.

Higher Volume

Maximum fill: **1325mL**

Volume markings reach 1000mL.

Shatter-Resistant

Lightweight plastic eliminates glass shard risk.



SV1000-SSP is the

Best of Both Worlds

between glass and plastic bottles.

Aspirate, Sample, and Remove Fluid Your Way.

REF #	ITEM	CASE QTY
SV1000-SSP	Spike-able Sterile Pathway 1000+mL Safety-Vac™ Bottle	12



NOT MADE WITH NATURAL RUBBER LATEX



NOT MADE WITH DEHP



NOT MADE WITH PVC



Ordering Information

Safety-Vac™ Evacuated Drainage Bottles

REF #	ITEM	CASE QTY
SV1000-SSP	Spike-able Sterile Pathway 1000mL Evacuated Drainage Bottle	12
SV1000-S	Standard Sterile 1000mL Evacuated Drainage Bottle + Tubing	12
SV1000-NS	Standard Sterile Pathway 1000mL Evacuated Drainage Bottle	12
SV1000-NSWC	Valved Standard Sterile Path. 1000mL Evacuated Drainage Bottle	12

Safety-Vac™ Evacuated Drainage Bottles are intended for use to aspirate, remove, or sample body fluids.

Drainage Accessories

REF #	ITEM	CASE QTY
SVHFSP-S	High-Flow Sterile 47" inch. Drain Line w/ Optional Spike Attachment	50
SVCT47-S	Standard Sterile 47" inch. male-male Drain Line	50

Jigsaw Medical® Drainage Accessories intended use: To aspirate, remove, or sample body fluids using an external vacuum source.



For all order inquiries, pricing, or product questions, please contact your Jigsaw Medical® representative.

Tel (630) 587-0019 **Fax** (630) 524-9088

Email info@jigsawmed.com



Refer to Instructions for Use for relevant warnings, precautions, complications, and contraindications.



Jigsaw Medical®
1012 Airpark Dr. Suite 1
Sugar Grove, IL 60554
www.jigsawmed.com
(630) 587-0019